

SCHILLER
The Art of Diagnostics



CARDIOVIT AT-2 PLUS

Two devices in one

Two devices in one



SCHILLER was the first manufacturer to combine the two most common pulmonary function and cardiac tests in one device - the CARDIOVIT AT-2 plus with ECG and spirometry function.

The CARDIOVIT AT-2 plus is a 6-channel ECG recorder with all (12) ECG channels simultaneously processed to provide instant ECG recordings. The device distinguishes itself through its absolute reliability, easy to use and an excellent cost-benefit ratio or the best price / value.


SIMPLE HANDLING - DOUBLE BENEFIT

The AT-2 plus includes the following features:

- Two tests in one device pays off
- Both functions are controlled from the same device
- Low weight and compact dimensions
- Large A4 size printout from integrated quality thermal printer
- Built-in rechargeable battery for mains-independent use - 4hrs normal use or 300 printouts on one battery charge
- Clear LCD for ECG preview prior to printing
- Simple one key operation for main functions
- Automatic or manual recording modes
- Selectable printing formats
- ECG memory for easy and unlimited number of copies
- Interpretation program option (including measurements) for children and adults
- The alphanumeric, water resistant keyboard with function keys allows easy and quick patient data entries and clinical remarks

INTERNATIONALLY RENOWNED ECG TECHNOLOGY:


- 6 channel real-time printout at 5, 25 or 50 mm/s
- Automatic 6 /12 channel printout on one or two pages
- 12 leads Rhythm ECG
- Quick and reliable PC based interpretation and measurement program for adult and children
- No need to repeat ECG recordings since the quality of the waveforms is optimized by the SCHILLER smoothing and baseline filters
- Save time and paper by simply monitoring all 12 leads on the integrated 3 channel screen



CARDIOVIT AT-2 PLUS

Two devices in one

SCHILLER
A M E R I C A S



FULFILLING SPIROMETRY STANDARDS:

The Spirometry function of the SCHILLER AT-2 PLUS is automatically available when the SP-250 sensor is plugged into the RS-232 connector on the left hand side of the unit.

The SP-250 includes the following features:

- Small and light sensor
- Inspiratory and expiratory pulmonary function tests: FVC, SVC, MV, MVV as well as pre/post medication tests
- Visual support due to real-time graphics and display of the measured values on the screen
- SP-250 sensor has a disposable single patient mouthpiece, the risk of cross contamination is eliminated

CARDIOVIT AT-2 PLUS

The background of the entire page is a repeating pattern of stylized hearts and lungs in a light gray color. The hearts are simple outlines, and the lungs are also simple outlines showing the main lobes. The pattern is dense and covers the entire area.

SCHILLER

A M E R I C A S

Call us: **1-888-845-8455**

Visit us at **schillerservice.com**

10903 NW 33rd St,
Doral, FL 33172, USA